



CORAVIN™

WINE ACCESS NEEDLES



The Coravin System uses a thin, non-coring medical grade needle to seamlessly pass through the cork of a bottle to access the wine. When the needle is removed, the cork reseals and continues to protect the wine. Each Wine Access Needle is durably constructed and Teflon™ coated, allowing it to withstand hundreds of insertions through the cork.

Features

- Medical Grade: Minimizes disruption to the cork, effectively allowing the cork to reseal completely after the needle is removed. Coravin's Wine Access Needle design was inspired by non-coring biopsy needles, which are delicate enough to be used during surgical procedures.
- Teflon™ Coating allows for effortless insertion and removal of the needle through the cork when using the Coravin System.
- Textured Grip Surface allows for quick removal and replacement of different needles for the Coravin System.
- Three sizes for use in a variety of situations and bottles.

Which wine access needle is right for me?

SPECIFICATION	Faster Pour Needle	Standard Needle	Vintage Needle
Gauge	16	17	18
Length	3.6"		
Band color (grip)	Red	Black	Silver
Materials	Teflon™ coated stainless steel		
When to use?	If you want a faster pour, and are accessing bottles with corks that are in good condition.	This is the ideal needle for most bottles and situations. It is also the same needle that comes standard with the Coravin System.	If you are accessing an older bottle, or are concerned about the quality of the cork.
Pour Speed (per 5oz glass)	16 to 20 seconds (About 24% faster than the standard needle)	± 21-25 Seconds	30 to 33 seconds (About 40% slower than the standard needle)

Product Specifications

SKU	Description	UPC	Boxed dimensions	Boxed weight
800302	Needle 16 Long Accy (Faster)	858976004115	1" x 4.75" x 1"	.02 lbs
800002	Needle 17 STD Accy	858976004054	1" x 4.75" x 1"	.02 lbs
800502	Needle 18 Long Accy (Vintage)	858976004122	1" x 4.75" x 1"	.02 lbs